at medical CERAMIC COATING



CERAMIC COATING

ORDER INFORMATION

Please add "C" to the end of the article number. e.g. F-0010C

... improves instruments:

- cost reduction and longer product life through
 higher surface hardness up to 10 times harder than stainless steel
 lower friction improves wear and tear of scissor blades
 more resistance against corrosion
- reduced eye-fatigue through less glare due to light absorbing black surface
- easier cleaning and reprocessing due to lower surface energy and non-stick effect
- anti-allergic due to absence of Chromium and Nickel in surface layer
- ... established in cardiology and neurology
- ... consists of TiAIN, applied in a magnetron sputtering process
- ... is 100% quality controlled and comes with lifetime warranty